



**Australian Office**

Ground Floor, 582 Heidelberg Rd,  
Fairfield, 3078  
Victoria, Australia  
T + 61 3 9329 1167  
F + 61 3 8732 0310  
ABN 43 120 905 839

**USA Office**

400 Continental Blvd  
El Segundo, CA, 90245  
USA  
T + 1 858 736 9160  
F + 1 858 346 9300

**UK Office**

18 King William Street  
London, United Kingdom  
EC4R 9AN  
T + 44 203 960 7512

[www.axxin.com](http://www.axxin.com)  
Commercial in Confidence

**Thank you for completing this form. Please send your response to [support@axxin.com](mailto:support@axxin.com)**

Axxin Support Form	
<b>Company:</b>	
<b>Name:</b>	
<b>Position:</b>	
<b>Email:</b>	
Instrument Issue	Desktop Software Issue
<p>The Instrument Serial Number is located on a label on the back side of the instrument.</p> <p><b>Instrument Serial Number:</b></p> <p><b>Revision Number:</b></p>	<p><i>Check the Desktop Software version number, go to: <b>Desktop Application / About</b></i></p> <p><b>Desktop software version:</b></p>
<p>If you are able to access the device application screens. Navigate to: <b>Settings / About</b>. You will find the entire device Software, Firmware and Hardware information here.</p> <p><b>Instrument Application Version x.x.x.x:</b></p>	<p><i>What Windows Operating System are you using?</i></p> <p><i>To find your system specification go to: <b>Control Panel / System and Security / System Information</b>. Record the Windows edition and the System type, e.g. Windows 7, 32 bit</i></p>
<p>In your own words describe what the problem is. On which screen you were looking at when the error occurred and if any USB devices such as a USB Memory key or Printer Attached. How frequently are you seeing the issue? (e.g. <i>intermittent - every few tests</i>)? Does the test run ever finish or does it always crash before this?</p> <p><b>Description of Fault:</b></p>	
<p><b>Please supply the any additional files such as an instrument/desktop application log file, screenshots, videos and photos.</b></p>	

PAGE NO.	FILE NAME	REVISION DESCRIPTION	TEMPLATE
1 of 1	D003976v3.0 Form, Support Form, Axxin Instrument Platform.docx	Added revision footer. Updated UK and US address	D002143v9.3